



The Naída CI Q Series Overview

November 2015



Advanced Bionics




What you will learn today

- The Naída CI Q70 and Q90 Processor Basics
- Indicators for System Checks and Troubleshooting
- Connectivity
 - Roger
 - ComPilot
 - Wireless Accessories





Naída CI Q Series





Naída CI Q Series

Child Friendly Features

- Light-emitting diode (LED) indicators
- Tamper proof to protect children from swallowing small parts
- Intellilink
- Ability to disable /limit volume and program settings to avoid accidental changes
- Flexible wearing options
- Waterproof with use of Aqua Case and Aqua Mic





What you will learn today

- **Processor Basics**
- Indicators for System Checks and Troubleshooting
- Care and Maintenance
- Connectivity
- Advanced Technology





Processor Basics

Program Button

- Stores up to five (5) programs
- Push the program button down briefly and then release the button to change programs
- The number of GREEN blinks displayed on the LED indicates which program is in use
- The programs switch in a chronological order; once it reaches the final program, it will return back to the first program
- The program button can be disabled by the child's audiologist

Program Button





Processor Basics

Did you know??

You do not need to change your student's programs for the different listening environments they encounter at school


Autosound OS™ processing **automatically** adapts to the changing listening environments at school, allowing children to hear every sound from whispers in the library to friends talking in the lunch room



Processor Basics

Volume Control

- To increase the volume, briefly tap the top of the volume control. To decrease the volume briefly tap the bottom of the volume control.
- The audiologist has the ability to restrict or disable the volume range





Processor Basics

Did you know??

The Naída CI will always default to program 1 at the child's standard volume setting when it is powered on.

If P1 with standard volume is the child's everyday program no setting changes are necessary to the Naída CI after the battery is attached!



Processor Basics

Microphones


- The Naída CI has four microphones
- The audiologist can set each program to utilize the microphones in different configurations

Headpiece Mic

Rear

Front

T- Mic 2



Check with the child's family or audiologist to determine which microphones are active in each program



Processor Basics- Batteries



PowerCel
110 mini



PowerCel
170 mini



PowerCel
170



PowerCel
230



Zinc Air

New



Processor Basics



Battery Type	Average (hours)
2 Zinc-Air Cartridge	31
110 PowerCel Mini	12
170 PowerCel Mini	18
170 PowerCel with FM	12
170 PowerCel without FM	18
230 PowerCel	25
AAA PowerPak	129



What you will learn today

- The Naída CI Q70 and Q90 Processor Basics
- **Indicators for System Checks and Troubleshooting**
- Care and Maintenance
- Connectivity
- Advanced Technology





Indicators for System Checks and Troubleshooting

Tri-Color LED

FEATURE	COLOR
Battery Status	Orange
Microphone Status	Green
Program Position	Green
CI Status	Red





Indicators for System Checks and Troubleshooting

ORANGE BEHAVIOR	BATTERY STATUS
Blinks at start up	<ul style="list-style-type: none">• 4 quick blinks battery is fully charged• 2-3 indicates that the battery is sufficiently charged to power the Naída CI• 1 quick blink indicates the battery is nearly depleted• No blinking indicates depleted battery
Solid	The battery is almost depleted
Blinks twice every 3 seconds	The battery is almost depleted and cannot support stimulation
Fades Out	The Naída CI is entering Standby Mode



Indicators for System Checks and Troubleshooting

GREEN BEHAVIOR	MICROPHONE/PROGRAM POSITION
Flickers in response to loud input	The sound processor and microphone are responding to sound
Blinks at start-up after battery status and upon program change	<ul style="list-style-type: none">• 1 blink indicates program one• 2 blinks indicates program two• 3 blinks indicates program three• 4 blinks indicates program four• 5 blinks indicates program five
Solid	A processor that is not yet programmed
Blinks 4 Times	Response to the AB myPilot request to "Find Paired Devices". The Left paired device will identify itself with this LED pattern.



Indicators for System Checks and Troubleshooting

RED BEHAVIOR	CI STATUS
Blinks once per second	Loss of lock with the implant.
Blinks rapidly (more than once per second)	Intellink™ is enabled and the Naída CI is connected to the wrong implant
Solid	Sound Processor Error Condition. Fully remove and re-insert the battery to reset processor
Blinks 5 Times	Response to AB myPilot request to “Find Paired Devices”. The Right paired device will identify itself with this LED pattern



Indicators for System Checks and Troubleshooting

Internal Alarms

ALARM	INDICATION
Beeps upon program change	1 beep indicates program 1 2 beeps indicates program 2 3 beeps indicates program 3 4 beeps indicates program 4 5 beeps indicates program 5
Short beep upon increase or decrease of volume	Beeps once per press of the volume control either up or down (a double beep will be heard when the following settings are reached: top of volume range; baseline setting; bottom of volume range)
Long beep once every 15 minutes	Low Battery



Indicators for System Checks and Troubleshooting

Intellilink

- Safety Feature
- Ensures the Naída CI only stimulates the correct internal device
- Prevents accidental stimulation from a wrong ear processor for a bilateral user or processor that belongs to another child
- The Naída CI LED will blink **RED** rapidly (more than once per second) to indicate it is connected to the wrong internal device



Indicators for System Checks and Troubleshooting

Naída CI Listening Check

- Verifies the Naída CI sound sources are working properly when children are unable to provide this feedback themselves
- Any sound source that is programmed for use with the child's Naída CI can be verified
- Sound Sources can include,
 - 4 Naída CI microphones
 - Wireless (via ComPilot)
 - FM/Roger
 - TCoil





Indicators for Status Check and Troubleshooting

Data logging

- Allows the audiologist to see data about how the Naída CI is being used on a daily basis
- The Audiologist can use this information to make programming changes to help the child maximize use of their device
- Can obtain information about
 - What programs and features are being used
 - What volume setting is being used
 - The different listening environments the child encounters throughout the day



AB MAKES IT SIMPLE FOR SCHOOLS!

Did you know??

AB recommends you obtain the following information from the child's audiologist to effectively manage your student's equipment in the classroom:

1. What program should be used for everyday listening?
2. What program should be used for FM/Roger or T-Coil use?



Troubleshooting

Basic Steps

When a child is not responding as expected:

1. Verify the UHP is in place and on the child's head
2. Remove the Naída CI and UHP from the child. Visually inspect the equipment and replace damaged parts
3. Verify the battery is charged
4. Reset the Naída to Program 1
5. Place the Naída CI and UHP back on the child and use the LED's and/or internal alarms to determine what the problem may be
6. Perform a listening check



Troubleshooting

Visually inspect the child's equipment and systematically replace damaged parts

- Check the Naída CI for damage
- Inspect the headpiece cable for any twisting, fraying, or breakage
- Verify the cable clicks or snaps into place when connected to the UHP
- Verify there is no visible damage to the UHP
- Inspect cable ports and jacks for debris. Clean with compressed air if needed
- Inspect the T-Mic 2 for any twisting, fraying, or breakage



Troubleshooting

Verify the battery is charged

- Remove the battery and then reattach it to the processor
- The LED (located in the middle of the volume control) will flash ORANGE to indicate battery status
- 3-4 ORANGE blinks indicates the battery is sufficiently charged to power the processor
- Replace with a fully charged battery if needed



Note: Zinc Air batteries will not provide LED battery status information upon start up. Only Power Cels and the AAA PowerPak will provide LED battery status information.



Troubleshooting

4. Re-set the child's processor to Program 1

- Remove and reattach the battery
- The Naída CI will always default to Program 1 with the volume set at the child's standard setting when the battery is removed and reattached
- If Program 1 is not the child's standard program, push the Program Button down briefly and then release the button until you reach the child's standard program

Note: If you have removed and reattached the battery as instructed in step 3 then the Naída CI has been re-set to Program 1.



Troubleshooting Situations

Troubleshooting Situations

Instructions for troubleshooting are available in the TFS Troubleshooting Guide

- No Sound and/or Red LED is flashing once per second
- No sound and/or there is a solid RED LED
- Sound has static, is muffled, or distorted
- Debris on battery contacts
- Naída CI does not power up
- Green LED does not flash in response to loud sounds
- No orange blinks or only one orange blink during battery LED sequence





Troubleshooting Situations

No Sound and/or RED LED is flashing once per second

1. Verify the UHP is positioned properly on the head
2. Remove any hats, scarfs, or other items that may be covering the microphone
3. Visually inspect the cable for any damage and verify it is firmly attached to the UHP and Naída CI
4. Remove and reattach the battery
5. Replace the UHP cable
6. Replace the UHP
7. Perform a listening check and replace any malfunctioning equipment



Troubleshooting Situations

No Sound and/or there is a solid RED LED

1. Remove and reattach the battery
2. Verify a charged PowerCel or two fully charged high power cochlear implant plus 675 Zinc-Air batteries are in place
3. Replace the UHP cable
4. Replace the UHP
5. Try a different program



Troubleshooting Situations

Static, Muffled, or Distorted Sounds

1. Remove any hats, scarfs etc. that may be covering the microphone
2. Verify the UHP is positioned properly on the head
3. Remove and reattach the battery
4. Visually inspect the cable for any damage and verify it is firmly attached to the UHP and Naída CI
5. Verify the Naída CI is set to the proper program and volume setting
6. If available, use an AB myPilot do a device status check and confirm program and volume settings
7. Try a different program
8. Visually inspect microphones for signs of debris and wear
9. Perform a listening check
10. Replace the UHP cable
11. Replace the UHP
12. Replace the T-Mic 2
13. Clean the battery contacts with compressed air



Troubleshooting Situations

Debris on Battery Contacts

If the battery contacts have rust or debris

1. Clean the contacts with compressed air
2. Place PowerCel batteries in the Zephyr Dry & Store[®] overnight
3. PowerCel batteries should be placed in the Zephyr Dry & Store[®] when not being charged





Troubleshooting Situations

Naída CI Does Not Power Up

1. Remove and reattach the battery
2. Verify the PowerCel is properly attached or the Zinc-Air batteries are inserted correctly
3. Verify you are using a fully charged PowerCel or two fully charged high power cochlear implant plus 675 Zinc-Air batteries

Note: Zinc-Air batteries will not provide LED battery status information upon startup.



Troubleshooting Situations

Green LED Does Not Flash in Response to Loud Sounds

1. The Naída CI must have LED's enabled by the audiologist
2. Remove and reattach the battery
3. Verify the PowerCel is properly attached or the Zinc-Air batteries are inserted correctly
4. Verify you are using a fully charged PowerCel or two fully charged high power cochlear implant plus 675 Zinc-Air batteries
5. Verify the processor is set to the proper program and volume setting
6. Try a different program
7. Perform a listening check
8. Replace the T-Mic 2
9. Clean the battery contacts with compressed air





Troubleshooting Situations

No Orange Blinks or Only One Orange Blink is Observed During Battery Status Check

1. Zinc-Air batteries will not provide LED battery status information upon start up. Only PowerCel batteries and the AAA PowerPak will provide LED battery status information
2. Remove and reattach the battery
3. Verify the PowerCel is properly attached to the Naída CI
4. Verify you are using a fully charged PowerCel
5. Clean the battery contacts with compressed air



What you will learn today

- The Naída CI Q70 and Q90 Processor Basics
- Steps for Everyday Use
- Indicators for System Checks and Troubleshooting
- Care and Maintenance
- **Connectivity**
- Advanced Technology





Connectivity

Connectivity Options

- T-Mic 2
- ComPilot
- FM/Roger
- T-Coil
- DECT Phone
- Easy Call





Connectivity

Natural Connectivity with the T-Mic 2

- Natural placement for headset, ear buds, or cell phones
- Access to commercially available audio devices





Connectivity

Wireless FM/Roger

Phonak ComPilot



Can be used with
Phonak's MLxi or Roger
X receivers

170 mAH battery



Can be used with Phonak's
Roger 17 receiver



ComPilot

- High quality audio streaming
 - Uses Phonak HIBAN proprietary technology
 - No t-coil antenna
- Wireless connectivity
 - MP3 players
 - TVs
 - Bluetooth Devices such as cell phones, laptop computers, automobiles
 - FM systems
- Remote Control functionality





What can ComPilot do?

Stream
Audio



Remote
Control



Stream
Roger/FM



Stream Bluetooth
Cell Phone or
Landline Phone



* Compatible with Phonak MLxi



ComPilot for Bimodal Streaming





Logical Button Arrangement

MAIN

- Start/ Pause Streaming
- Switch programs

VOLUME

- If programmed as a remote, it controls the CI volume only.
- If not programmed as a remote, it controls the streaming volume only.

CONNECT

- Used with + button, holding for 2 seconds to begin Bluetooth pairing

Home Orientation

- Short press toggles through programs
- Long press returns to program 1 with default volume settings



- Europort to connect MLXi



What is Roger?



Roger is the new digital standard that bridges the understanding gap, in noise and over distance, by wirelessly transmitting the speaker's voice directly to the listener.



How does Roger work?

- Hearing Assistive Technology
 - Superior to traditional and Dynamic FM
 - For hearing aids, cochlear implants and/or soundfield amplification systems
- Adaptive Digital Wireless Technology*
 - Utilizes digitized audio “packets”
 - Sent on different channels (2.4000-2.44835 GHz)
 - Hop between channels
 - Only use free channels
 - Two-way communication between the microphone and receiver to ensure privacy for the user



Benefits?

- No interference
- No noticeable transmission delay (below 25 ms)
- No limit to the number of systems which can be used in one building
- No need to select specific channels for each classroom




Research on Roger and CIs

- Jace Wolfe, Hearts for Hearing
- 37 participants
- 16 AB recipients using Harmony
- Roger systems in combination with CIs resulted in significant improvements in speech recognition at high noise levels (70, 80 dB) over traditional and Dynamic FM technologies






Research on Roger and CIs



- No significant difference in quiet
- No significant difference between Roger, MLxi, or MLxS in low level noise



Research on Roger and CIs*



- At higher noise levels (70, 75, and 80 dB) Roger provided significantly better speech recognition in noise when compared to traditional FM
- Roger provided better performance than Dynamic FM



Roger Receivers for Naída CI



Roger myLink



Roger 17



Roger X*

* Roger X is compatible with Naída CI via Phonak ComPilot



Roger 17

- Design integrated for Naída CI Q70
- Connects to the bottom of the 170 PowerCel
- 1 gram light, 8 mm short




170 PowerCel





Roger 17 (02 & 03)

Compatibility between microphones and receivers





Roger Pen

- Adaptive wireless transmission
 - Automatic microphone mode selection
- Bluetooth for cell phone use
- TV connectivity
- Audio input for listening to multimedia
- Used with Roger Clip-On Mics and Roger Pens in a microphone network





Roger Pen






Roger Pen

The Roger Pen features an accelerometer

- An accelerometer measures accelerations in three dimensions (X, Y and Z).
- The accelerometer informs the Roger Pen continuously about its orientation with respect to the direction of gravity.
- It also tells the Roger Pen when it is accidentally dropped.





Roger Docking Station

Listening to multimedia via docking station

- Automatic detection of an audio signal (even when OFF)
 - Powers down if no audio signal for 45 seconds
 - No power down even during soft music





Clip-On Mic for Roger (02 & 03)

- Continuous analysis of surrounding speech and noise
- Automatic configuration for beamforming, noise suppression and gain settings automatically
- No Bluetooth connectivity
- Audio cable to connect to TV and other audio devices
- Allows for multiple microphone use to create a MultiTalker Network





Clip-On Mic for Roger (02 & 03)

Indicator light

Microphones

Charging and Audio input (micro-USB)



On / Off and Mute



Connect





Roger Resources

Fitting Guide

roger

Roger and Cochlear Implants

This guide provides detailed information on how Roger should be used with the most popular cochlear implants to achieve the best possible performance.

A recent study by Dr. Jace Wolfe of the University of Oklahoma, Oklahoma City, revealed that combination with cochlear implants resulted in significant improvements in speech recognition over fixed gain FM and Dynamic FM technologies (see www.phonakpro.com).

Set-up

The table below shows what Roger receiver and/or adapter is required to use Roger with a given product name.

Advanced Bionics	Advanced Bionics	Product name
Naída CI Q70	Harmony® / Auro®	Cochlear Nucleus 6 (Cp6) / Nucleus 5 (Cp5)
Roger 17™ or ComPilot™ + Roger X	iConnect™ + Roger X ²	Roger 14 or adapter 2

Roger X is not compatible with the Advanced Bionics (AB) Neptune™ sound processor.

¹ Jace Wolfe (2015) Evaluation of speech recognition of cochlear implant recipients using a listening system. Accepted by the Journal of the American Academy of Audiology.

² PowerMax™ 170 battery required.

³ Roger X with SR-1350N3500 only.

⁴ Build Standard C or higher.



Device	Receivers	Microphones
Naída CI via Direct Audio Input (DAI)	Roger 17	Roger inspiro Roger Pen Roger Clip-On Mic
Naída CI via Telecoil	Roger myLink	Roger inspiro Roger Pen Roger Clip-On Mic
Naída CI via ComPilot	Roger X	Roger inspiro Roger Pen Roger Clip-On Mic
Harmony via iConnect	Roger X	Roger inspiro Roger Pen Roger Clip-On Mic
Harmony via T-Coil	Roger myLink	Roger inspiro Roger Pen Roger Clip-On Mic
Neptune via T-Comm™	Roger myLink	Roger inspiro Roger Pen Roger Clip-On Mic

* Roger DynaMic is a pass-around microphone that works alongside Roger inspiro in a classroom's MultiTalker Network.

Fitting Remote Microphone Technology for AB Recipients

Naída CI Q70

Audio Mixing	Additional Program Settings
Need Aux input (50/50, 30/70, Aux only)	- Internal Telecoil disabled - Mic Mode set to Omnidirectional
No need for Aux input	- Internal Telecoil disabled - Mic Mode set to Omnidirectional
Need Aux input (50/50, 30/70, Aux only)	- Set ComPilot percentage in Processor Pane prior to download (25%, 50%, 75%, 100%) - Internal Telecoil disabled - Mic Mode set to Omnidirectional

70 via the PowerCell 170 battery.

17, identifiable by the numbers 02 or 03 on the receiver.

1 with all Roger microphones (Roger inspiro, Roger Pen, and Roger Clip-On Mic only).

with Roger Pen and Roger Clip-On Mic only.

Use of Roger Pen and Roger Clip-On Mic.

Roger inspiro	Roger Pen	Roger Clip-On Mic
Yes	No	No
No	Yes	No
Yes	Yes	Yes

Audio Mixing

Audio Mixing	Additional Program Settings
Need Aux input (50/50, 30/70, Aux only)	- Aux Source set to Auto-Detect or FM only - T-comm switch in middle or bottom position
Aux input (50/50, 30/70, Aux only)	- Internal Telecoil disabled
Aux input (50/50, 30/70, Aux only)	- Internal Telecoil disabled

TOOLS for SCHOOLS™



The POWER of TWO



The Power of Two
Working Together as One